

Lunit INSIGHT CXR

AI Solution for Chest X-ray





Lunit INSIGHT CXR

Product Information

10 Product Overview

14 User Benefits

18 Clinical Validation

22 Clinical Cases

UK • Netherlands

China •
South Korea

• United States

• Mexico


• Brazil


200⁺
Reference Sites

80⁺
Countries Worldwide

7M⁺
*Images Analyzed
(for clinical and research use)*

Never miss
a finding.

Lunit INSIGHT CXR LOW  HIGH



CASE REPORT

Abnormality Score 99%

Atl	Atelectasis	98%	Multiple Lesions
PEf	Pleural effusion	85%	Right zone
Calc	Calcification	Low	-
Cm	Cardiomegaly	Low	-
Csn	Consolidation	Low	-
			⋮

What does Lunit INSIGHT CXR analyze on chest x-ray images?

Lunit INSIGHT CXR detects 10 abnormal radiologic findings with 97-99% accuracy and supports tuberculosis screening on chest x-ray images.

10 *Abnormal Radiologic Findings*

97-99% *Accuracy*



Lunit INSIGHT CXR generates

Developed by using Lunit's cutting-edge deep learning technology


Detected Location

The location information of detected lesions in the form of heatmaps and/or contour maps



Abnormality Score

An abnormality score reflecting the AI's calculation of the actual presence of the detected lesion



AI Report

An AI case report summarizing the overall analysis result, narrowed down to each finding

CASE REPORT

Abnormality Score 99%			
Atel	Atelectasis	99%	Multiple Lesions
PEf	Pleural effusion	88%	Right zone
Calc	Calcification	Low	-
Cm	Cardiomegaly	Low	-
Csn	Consolidation	Low	-
			⋮

What are the major benefits of using it?

Fast triage of normal cases.

Efficient reading via exam prioritization.

Improved reading performance.

01

Fast triage of normal cases

According to the abnormality scores generated by AI, radiologists can triage normal cases quickly and focus on reading abnormal cases where lesions might exist.

Health check-up centers

Imaging clinics

Teleradiology centers

02

Efficient reading via exam prioritization

In reference to the abnormality scores on the worklist, radiologists can prioritize exams in their reading order, resulting in a 13% reduction in reading time, and a 33% in reduction time for normal cases.¹

Radiology departments

Imaging clinics

Teleradiology centers

13%
*reduction
in reading time*

33%
*in reduction time for
normal cases*

03

Improved reading performance

Non-radiology physicians, general radiologists, and thoracic radiologists can improve their diagnostic accuracy for major chest abnormalities such as malignant pulmonary nodules, pneumothorax, pneumonia, and active pulmonary tuberculosis.^{2 3 4 5 6 7 8}

Clinical departments prior to surgeries and procedures

(Emergency, Pulmonary, Oncology, Cardiothoracic, Surgery, etc.)

Community hospitals and clinics

Reduced overlooked lung cancers.

Streamlined ED workflow.

COVID-19 patient triaging and monitoring.

04

Reduced overlooked lung cancers

The AI-aided, automatic detection of small and subtle pulmonary nodules overlapped in the hilar shadow, ribs, heart, and diaphragm, enables radiologists to reduce overlooked lung cancer cases, especially during regular check-ups.⁹

Health check-up centers

Community hospitals and clinics

05


Streamlined ED workflow

With the assistance of AI, radiology residents can improve their diagnostic performance¹⁰ and reduce their reading time¹¹, which ultimately accelerates the decision-making process and treatment in the ED.

Emergency departments

<Reading Time for Detecting Acute Respiratory Infection>

(Image N=387)



06

COVID-19 patient triaging and monitoring

AI-aided chest radiograph interpretation can help medical professionals detect COVID-19 infected pneumonia quickly¹² and accurately¹³, enabling prompt isolation and timely treatment.

COVID-19 screening centers

Emergency departments

- [Click to Watch](#)
- [Thammasat University Hospital, Thailand](#) [↗](#)
- [Songklanagarind Hospital, Thailand](#) [↗](#)

What do the medical journals say?


Below are highlights from the studies published in peer-reviewed journals that validate the performance of Lunit INSIGHT CXR and its clinical value in chest radiograph interpretation.



Accurate and efficient diagnosis boosted with AI


HIGHLIGHT 1

Improved reading performance of non-radiology physicians, general radiologists, and thoracic radiologists.¹⁴




HIGHLIGHT 2

Better detection of early-stage overlooked lung cancer without increasing false positive cases.¹⁵




HIGHLIGHT 3

Reduce the overall reading time.¹⁶



HIGHLIGHT 4


AI-aided chest radiograph interpretation takes less time than PCR tests in triaging COVID-19 suspected patients.¹⁷



<Turnaround Time Between Radiology Report and rRT-PCR Result>


How does AI work in COVID-19 settings?

Go to Video [Testing COVID-19](#)



What do the medical journals say about AI-powered chest x-ray interpretation?


Go to Video [→](#)



Five-year survival rate when detected early by AI

Lunit INSIGHT successfully analyzed the chest x-ray image of 54-year-old male patient, detecting lung cancer that had been missed 3 years ago.

73%
When Detected Early By AI (stage 1-2)



18%
When Missed (stage 3-4)



Reference : AJCC 8th Edition

2013



MISSED




AI Score 16.7%
DETECTED

2014



MISSED




AI Score 43.1%
DETECTED

2016




LUNG CANCER
DIAGNOSED




AI Score 90.7%
DETECTED

Sample Cases




Small nodule detected

63%
Abnormality Score



Multiple lesions including atelectasis, pleural effusion

98%
Abnormality Score



COVID-19 infected pneumonia detected

77%
Abnormality Score

Try Lunit AI Solution

Visit insight.lunit.io and upload a DICOM file you have. You will get the AI result within seconds.

Go to page [→](#)



Reference

☰ User Benefits

- ¹ Ju Gang Nam, Minchul Kim, et al. Development and validation of a deep learning algorithm detecting 10 common abnormalities on chest radiographs. *European Respiratory Journal*. 2020
- ² Ju Gang Nam, Sunggyun Park, et al. Development and Validation of Deep Learning–based Automatic Detection Algorithm for Malignant Pulmonary Nodules on Chest Radiographs. *Radiology*. 2018
- ³ Eui Jin Hwang, Sunggyun Park, et al. Development and Validation of a Deep Learning-based Automatic Detection Algorithm for Active Pulmonary Tuberculosis on Chest Radiographs. *Clinical Infectious Diseases*. 2018
- ⁴ Eui Jin Hwang, Sunggyun Park, Kwang-Nam Jin, et al. Development and Validation of a Deep Learning–Based Automated Detection Algorithm for Major Thoracic Diseases on Chest Radiographs. *JAMA Network Open*. 2019
- ⁵ Jong Hyuk Lee, Sunggyun Park, et al. Deep learning–based automated detection algorithm for active pulmonary tuberculosis on chest radiographs: diagnostic performance in systematic screening of asymptomatic individuals. *European Radiology*. 2020
- ⁶ Eui Jin Hwang, Jung Hee Hong, et al. Deep learning algorithm for surveillance of pneumothorax after lung biopsy: a multicenter diagnostic cohort study. *European Radiology*. 2020
- ⁷ Jong Hyuk Lee, Hye Young Sun, et al. Performance of a Deep Learning Algorithm Compared with Radiologic Interpretation for Lung Cancer Detection on Chest Radiographs in a Health Screening Population. *Radiology*. 2020
- ⁸ Hyunsuk Yoo, Ki Hwan Kim, et al. Validation of a Deep Learning Algorithm for the Detection of Malignant Pulmonary Nodules in Chest Radiographs *JAMA Network Open*. 2020
- ⁹ Sowon Jang, Hwayoung Song, et al. Deep Learning–based Automatic Detection Algorithm for Reducing Overlooked Lung Cancers on Chest Radiographs. *Radiology*. 2020
- ¹⁰ Eui Jin Hwang, Ju Gang Nam, et al. Deep Learning for Chest Radiograph Diagnosis in the Emergency Department. *Radiology*. 2019
- ¹¹ Jae Hyun Kim, Jin Young Kim, et al. Clinical Validation of a Deep Learning Algorithm for Detection of Pneumonia on Chest Radiographs in Emergency Department Patients with Acute Febrile Respiratory Illness. *Journal of Clinical Medicine*. 2020
- ¹² Eui Jin Hwang, Hyungjin Kim, et al. Implementation of a Deep Learning-Based Computer-Aided Detection System for the Interpretation of Chest Radiographs in Patients Suspected for COVID-19. *Korean Journal of Radiology*. 2020
- ¹³ Se Bum Jang, Suk Hee Lee, et al. Deep-learning algorithms for the interpretation of chest radiographs to aid in the triage of COVID-19 patients: A multicenter retrospective study. *PLOS ONE*. 2020

☰ Clinical Validation

- ¹⁴ Eui Jin Hwang, Sunggyun Park, Kwang-Nam Jin, et al. Development and Validation of a Deep Learning–Based Automated Detection Algorithm for Major Thoracic Diseases on Chest Radiographs. *JAMA Network Open*. 2019
- ¹⁵ Sowon Jang, Hwayoung Song, et al. Deep Learning–based Automatic Detection Algorithm for Reducing Overlooked Lung Cancers on Chest Radiographs. *Radiology*. 2020
- ¹⁶ Ju Gang Nam, Minchul Kim, et al. Development and validation of a deep learning algorithm detecting 10 common abnormalities on chest radiograph, *European Respiratory Journal*, 2020
- ¹⁷ Eui Jin Hwang, Hyungjin Kim, et al. Implementation of a Deep Learning-Based Computer-Aided Detection System for the Interpretation of Chest Radiographs in Patients Suspected for COVID-19. *Korean Journal of Radiology*. 2020

Research Partners

We welcome research partnerships and other collaboration with medical institutions, healthcare providers, and companies interested in implementing our software products. Currently, we have over 30 worldwide research partners throughout USA, Canada, UK, Sweden, China, and Korea.

Distribution Partners

FUJIFILM



AGFA *Agfa*

SECTRA

PHILIPS



This document is for use by healthcare professionals only. The radiologist should always rely on his own clinical and professional opinion when deciding whether to use a certain product to diagnose or treat a patient.

Not all Lunit products may be available in all markets, as availability is based on the medical and/or regulatory practices of each market. Contact your Lunit representative if you have questions about the availability of the Lunit products in your area.

Lunit Inc. owns, uses or enforces the following trademarks or service marks: Lunit, Lunit INSIGHT CXR. The products described are CE-marked in accordance with applicable EU standards and directives. The product fulfils the requirements for CE marking.

Lunit Inc. 15 Floor, 27 Teheran-ro 2-gil, Gangnam-gu, Seoul, 06241, Republic of Korea

FDA clearance expected in 2021
Document number: CXR-BR-TB-EN_Ver1

Contact us

General - contact@lunit.io
Investment - ir@lunit.io
Media - media@lunit.io
Partnership - partner@lunit.io

Office

Seoul HQ
Boston, USA
Amsterdam, The Netherlands
Shanghai, China

AI Vision, Earlier Action

Contact Us

Please feel free to email us about any inquiries or questions.

E-mail : contact@lunit.io

Website : www.lunit.io