




**APERIO® Hybrid<sup>17|21</sup>**  
Thrombectomy Device



For vessel diameters from 1.0 – 5.5 mm  
Effective hybrid cell design  
Visible throughout the entire length

# ORDERING INFORMATION | APERIO® Hybrid<sup>17|21</sup>

Labelled APERIO® Hybrid <sup>17 21</sup> Dimensions (mm)	Reference Number	Device Diameter (mm)	Device Length* (mm)	Recommended Vessel Diameter (mm)	Required / Recommended Microcatheters for Delivery (Inch)
2.5 × 16	01-000713	2.5	16	1.0 – 2.0	0.0165 – 0.021 NeuroSlider® 17 DLC NeuroSlider® 21 DLC
2.5 × 28	01-000710	2.5	28	1.0 – 2.0	
3.5 × 28	01-000711	3.5	28	1.5 – 3.0	
4.5 × 30	01-000712	4.5	30	2.0 – 4.0	
4.5 × 40	01-000715	4.5	40	2.0 – 4.0	0.021 – 0.027 NeuroSlider® 21 DLC NeuroSlider® 27 (DLC)
4.5 × 50	01-000716	4.5	50	2.0 – 4.0	
6.0 × 40	01-000717	6.0	40	3.5 – 5.5	
6.0 × 50	01-000718	6.0	50	3.5 – 5.5	

\* Average length within intended vessel diameter

All changes or modifications, may they be technical or other, or changes in the availability of products are expressively reserved.

