



**CARESTO® heal**  
Stent



For treatment of carotid artery stenosis

Antithrombogenic coating \*


Different device configurations for tailored treatment options

Perfect flexibility and high plaque coverage


\* The antithrombogenic properties were demonstrated in non-clinical tests.

# ORDERING INFORMATION | CARESTO® heal

| Labelled CARESTO® heal Dimensions (mm) | Reference Number | Unconstrained Device Diameter (mm) | Unconstrained Device Length (mm) | Recommended Vessel Diameter (mm) | Required/ Recommended Catheters for Delivery |
|--|------------------|------------------------------------|----------------------------------|----------------------------------|--|
| 6.0 × 30                               | 01-005301        | 6.3                                | 23                               | 4.0 – 6.0                        | NeuroSlider® 52 DLC (ID: 0.052")             |
| 6.0 × 40                               | 01-005302        | 6.3                                | 30                               | 4.0 – 6.0                        |  |
| 6.0 × 50                               | 01-005303        | 6.3                                | 37                               | 4.0 – 6.0                        |  |
| 8.0 × 30                               | 01-005311        | 8.3                                | 24                               | 6.0 – 8.0                        |  |
| 8.0 × 40                               | 01-005312        | 8.3                                | 32                               | 6.0 – 8.0                        |  |
| 8.0 × 50                               | 01-005313        | 8.3                                | 40                               | 6.0 – 8.0                        |  |
| 8.0 × 60                               | 01-005314        | 8.3                                | 47                               | 6.0 – 8.0                        |  |
| 10.0 × 30                              | 01-005321        | 10.3                               | 25                               | 8.0 – 10.0                       |  |
| 10.0 × 40                              | 01-005322        | 10.3                               | 33                               | 8.0 – 10.0                       |  |
| 10.0 × 50                              | 01-005323        | 10.3                               | 41                               | 8.0 – 10.0                       |  |
| 10.0 × 60                              | 01-005324        | 10.3                               | 49                               | 8.0 – 10.0                       |  |



| Labelled CARESTO® heal Dimensions (mm) | Reference Number | Unconstrained Device Diameter dist. / prox. (mm) | Unconstrained Device Length (mm) | Recommended Vessel Diameter dist. / prox. (mm) | Required/ Recommended Catheters for Delivery |
|--|------------------|--|----------------------------------|--|--|
| 6.0 / 8.0 × 29                         | 01-005351        | 6.3 / 8.3  | 29                               | 4.0 – 6.0 / 6.0 – 8.0                          | NeuroSlider® 52 DLC (ID: 0.052")             |
| 6.0 / 8.0 × 35                         | 01-005352        | 6.3 / 8.3  | 35                               | 4.0 – 6.0 / 6.0 – 8.0                          |  |
| 6.0 / 8.0 × 41                         | 01-005353        | 6.3 / 8.3  | 41                               | 4.0 – 6.0 / 6.0 – 8.0                          |  |
| 8.0 / 10.0 × 29                        | 01-005361        | 8.3 / 10.3                                       | 29                               | 6.0 – 8.0 / 8.0 – 10.0                         |  |
| 8.0 / 10.0 × 36                        | 01-005362        | 8.3 / 10.3                                       | 36                               | 6.0 – 8.0 / 8.0 – 10.0                         |  |
| 8.0 / 10.0 × 43                        | 01-005363        | 8.3 / 10.3                                       | 43                               | 6.0 – 8.0 / 8.0 – 10.0                         |  |



All changes or modifications, may they be technical or other, or changes in the availability of products are expressly reserved.  
Not available for sale in the United States.

CE 0297